



WriteSite® Plus Surgical Skin Marker

Intended Use

Intended to provide an accurate and reliable way to mark the incision site prior to surgery with Gentian Violet Ink formulated specifically for use with Chloraprep® Products.

Key Features

- Ink is three times darker than standard Gentian Violet ink after prep with Chloraprep® products
- Ink is also compatible with Duraprep® and Betadine® prep solutions
- Ink formulated to help reduce blurring, leaving the mark clear and visible for prep
- Allows for visibility of surgical markings to support correct site surgery and the Joint Commission's Universal Protocol
- Offered in a tip style designed specifically for precision marking
- Available in regular and small barrel sizes
- Also available as a multi-ink marker including black utility marker ink on the opposite end
- Ink requires only one minute of dry time before prep for optimal results
- Virtually eliminates the need for remarking after prep
- Ink is non-tattooing and non-irritating
- Intended for single-use and not made with natural rubber latex

WriteSite® Plus

Surgical Skin Marker

Formulated for use with ChloroPrep® Products.



WriteSite® Plus Standard Size Surgical Skin Marker

Technical Specifications

Part No.	Description	Sterility	Qty
2700	Marker Only	Sterile	50/Box
2701	Marker, Ruler & Label Set	Sterile	50/Box
2700BN	Marker Only	Non-Sterile	100/Box
2705BN	Marker Only	Non-Sterile	50/Jar

WriteSite® Plus Jr Size Surgical Skin Marker

Technical Specifications

Part No.	Description	Sterility	Qty
2715	Marker Only	Sterile	50/Box
2716	Marker, Ruler & Label Set	Sterile	50/Box
2715BN	Marker Only	Non-Sterile	100/Jar

WriteSite® Plus Multi-Ink Marker

Technical Specifications

Part No.	Description	Sterility	Qty
2720	Marker, Ruler & Label Set	Sterile	50/Box
2720BN	Marker Only	Non-Sterile	100/Box
2721BN	Marker Only	Non-Sterile	500/Box

New!

WriteSite® is a registered trademark of Aspen Surgical Products, Inc. ChloroPrep® is a registered trademark of CareFusion 2200, Inc. Duraprep® is a registered trademark of 3M Company. Betadine® is a registered trademark of Purdue Products L.P. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2021 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
M-SM-CAT-1000896 rev 13 23-Nov-2021 ENG – US

Aspen Surgical Products, Inc.
phone +1 616.698.7100 toll-free +1 888.364.7004
aspensurgical.com

