

# SD3

Ultrasonic Pocket Doppler



**EDAN**

# SD3 Ultrasonic Pocket Doppler

With excellent sensitivity, complete interchangeability and high durability, SD3 series Ultrasonic Pocket Dopplers are ideal for routine fetal heart rate detection by clinicians.

## User-Friendly Design

- Compact design with easy one-hand operation
- Bright OLED screen with real-time FHR display
- Hi-Fi speaker

## Interchangeable Waterproof Probe

- Plug-and-play
- 2 MHz probe with deep penetration for later pregnancy and large patients
- 3 MHz probe detects fetal heart beat as early as 9 weeks
- 4/5/8 MHz probe for blood flow detection
- Automatic power-off while putting the probe back to the rack (SD3 PLUS & SD3 PRO)

## Work Mode

- Real-time work mode for spot check
- Record & play with audio playback and numeric display

## Power Source

- Standard AA batteries for easy replacements
- Built-in Li-ion rechargeable battery ensures long working time
- Unique charger stand and wall mount design



# Configurations

Model	SD3 VASCULAR	SD3 LITE	SD3	SD3 PLUS	SD3 PRO
FHR Display	×	×	√	√	√
Charger Stand	×	×	×	○	○
Record & Play	×	×	×	×	√
AA Batteries	○	○	○	×	×
Li-ion Battery	×	×	×	√	√

√ Standard    × Not available    ○ Optional



SD3 VASCULAR



SD3 LITE



SD3



SD3 PLUS



SD3 PRO





**SD3**  
Ultrasonic Pocket Doppler



**Global Headquarters:**

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen  
518122 P.R. China | +86.755.26898326 | [www.edan.com](http://www.edan.com) | [info@edan.com](mailto:info@edan.com)

**U.S. and Canada inquiries:**

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126  
+1.858.750.3066 | [www.edandiagnostics.com](http://www.edandiagnostics.com) | [edan-info@edandiagnostics.com](mailto:edan-info@edandiagnostics.com)

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-OB & GYN  
SD3-V2.5-20240517