



Top Neuropole SU Needle

For a clear ultrasound image

New Bevel

for better punctures through the skin
and a safe nerve approach





• Ultrasound reflection with a regular plexus needle.

• Optimal ultrasound reflection with Top's Circle Area at the tip of the needle.

Top Neuropole SU Needle

Top combines over 30 years of experience developing plexus needles, with the latest ultrasound reflection technology. The echogenic Neuropole SU Needle, specifically designed for ultrasound-guided regional anesthesia procedures, provides an unequalled visibility.

The Circle Area at the tip of the needle generates enhanced ultrasound reflection. Even at a 60° angle or more, the Top Neuropole needle will provide you the clear image you need to steer the needle close to the nerves.

Insulated with Top's high grade medical Teflon-coating for smooth penetration, the Neuropole SU Needle will prove to be a great aid in optimizing your procedure, at a very attractive price level.

Ordering information

| TOP NEUROPOLE SU Needles, with electrical stimulation: | | | | |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|--------------|-------------------|----------------|
| <i>Reference number</i> | <i>Needle length</i> | <i>Gauge</i> | <i>Active tip</i> | <i>Package</i> |
| 627240 <i>(prev. 627181*)</i> | 30 mm | 23 G | 1 mm | 10 pcs/box |
| 627241 <i>(prev. 627182*)</i> | 50 mm | 23 G | | |
| 627242 <i>(prev. 627197*)</i> | 60 mm | 22 G | | |
| 627243 <i>(prev. 627198*)</i> | 100 mm | 21 G | | |
| 627244 <i>(prev. 627201*)</i> | 150 mm | 20 G | | |
| TOP NEUROPOLE RU Needles, without electrical stimulation: | | | | |
| 627245 <i>(prev. 627202*)</i> | 60 mm | 22 G | - | 10 pcs/box |
| 627246 <i>(prev. 627203*)</i> | 100 mm | 21 G | | |
| <small><i>Ad* The new Neuropole SU Sharp needles are replacing the original needles, with the "prev." reference numbers mentioned above. These older versions will however remain available on special request.</i></small> | | | | |

