



MEDICAL Z DLP TECHNOLOGY®

A UNIQUE TECHNOLOGY THAT COMBINES **COMPRESSION**
WITH **LYMPHATIC DRAINAGE**



Models registered and patent pending

DLP TECHNOLOGY® compression garments offer patients safety, efficiency and comfort in the post-operative period following surgery.

The compression they provide, reduces lymphatic stasis reducing postoperative pain, accelerates healing and creates optimal skin draping.

MEDICAL Z has developed an **MAJOR INNOVATION** in post-operative care for liposuction and abdominoplasty.

With a view to improving results,
Medical Z now offers
**A COMPRESSION GARMENT INCLUDING
AN INTEGRATED LYMPHATIC DRAINAGE SYSTEM CALLED :**
DLP TECHNOLOGY®



PRE-OPE	J+7 POST-OPE	J+15 POST-OPE	J+21 POST-OPE												
	<table border="1"><tr><td>TRADITIONAL COMPRESSION + DLP TECHNOLOGY®</td><td>TRADITIONAL COMPRESSION</td></tr><tr><td></td><td></td></tr></table>	TRADITIONAL COMPRESSION + DLP TECHNOLOGY®	TRADITIONAL COMPRESSION			<table border="1"><tr><td>TRADITIONAL COMPRESSION + DLP TECHNOLOGY®</td><td>TRADITIONAL COMPRESSION</td></tr><tr><td></td><td></td></tr></table>	TRADITIONAL COMPRESSION + DLP TECHNOLOGY®	TRADITIONAL COMPRESSION			<table border="1"><tr><td>TRADITIONAL COMPRESSION + DLP TECHNOLOGY®</td><td>TRADITIONAL COMPRESSION</td></tr><tr><td></td><td></td></tr></table>	TRADITIONAL COMPRESSION + DLP TECHNOLOGY®	TRADITIONAL COMPRESSION		
TRADITIONAL COMPRESSION + DLP TECHNOLOGY®	TRADITIONAL COMPRESSION														
TRADITIONAL COMPRESSION + DLP TECHNOLOGY®	TRADITIONAL COMPRESSION														
TRADITIONAL COMPRESSION + DLP TECHNOLOGY®	TRADITIONAL COMPRESSION														
SMOOTHER RESULTS WITH DLP TECHNOLOGY®.															



DO NOT WEAR IF WE HAVE A SILICONE ALLERGY

Models registered and patent pending

This **NEW CONCEPT** of innovative post-operative compression incorporates a lymphatic drainage system directly on the compression garment for an **ACCELERATED POSTOPERATIVE REHABILITATION.**

This significant technical improvement through drainage results in **REDUCED EDEMA AND ECCHYMOSIS,** for **FASTER OPTIMIZED RESULTS.**




EXCEPTIONAL RESULTS

From the study conducted by Professor Mojallal, head of the plastic, reconstructive and aesthetic surgery department, Croix-Rousse Hospital, CHU de Lyon, Hospices CIVILS DE LYON

Models registered and patent pending

The compression garment must not be worn directly on injured skin.



DRAIN BELT FULL TORSO DLP TECHNOLOGY® - DR/032			
SIZES	UNDER BREAST (cm) ③	WAIST CIRC. (cm) ④	HIP CIRC. (cm) ⑤
1	70 - 75	67 - 72	83 - 88
2	76 - 81	73 - 78	89 - 93
3	82 - 87	79 - 84	94 - 99
4	88 - 93	85 - 90	100 - 105
4+	94 - 99	91 - 96	106 - 112
5	100 - 105	97 - 102	113 - 119
6	106 - 115	103 - 112	120 - 130
7	116 - 121	113 - 118	131 - 141

DRAIN LIPO-PANTY DLP TECHNOLOGY® - DR/002				
SIZES	UNDER BREAST (cm) ③	WAIST CIRC. (cm) ④	HIP CIRC. (cm) ⑤	THIGH (cm) ⑥
1	71 - 75	63 - 67	83 - 88	44 - 49
2	76 - 80	68 - 72	89 - 93	50 - 55
3	81 - 85	73 - 77	94 - 99	55 - 59
4	86 - 90	78 - 82	100 - 105	59 - 62
4+	91 - 95	83 - 87	106 - 112	63 - 67
5	96 - 100	88 - 92	113 - 119	68 - 72
6	101 - 108	93 - 100	120 - 130	73 - 77
7	109 - 113	101 - 108	131 - 141	78 - 82

Lipo-panty® and combi : measures taken for patients between 1m60 and 1m75.



DO NOT WEAR IF WE HAVE A SILICONE ALLERGY

ZI N°1 - 14, rue Georges Cuvier - 37550 SAINT-AVERTIN
 Tél. : +33 (0)2 47 71 33 33 - Fax. : +33 (0)2 47 71 33 34
 Email : infofrance@medicalz.com

WWW.MEDICALZ.COM

Models registered and patent pending