



MEDICALZ®

CRÉATEUR FRANÇAIS DE VÊTEMENTS COMPRESSIONS DEPUIS 1981

AESTHETIC & PLASTIC SURGERY



WWW.MEDICALZ.COM

THE HIGHEST
QUALITY FOR YOUR
WELL-BEING





SUMMARY

INTRODUCING MEDICAL Z® 4

WOMEN'S GARMENTS:

Breast surgery:
Breast and body garments 7

Body contouring surgery:
Lipo-panty® 8
Lipo-panty®, Combi suit and Arm Sleeves 10
Combi® and Panty-lift® 11

Abdominal surgery:
Post-surgery belt 12

MEN'S GARMENTS 14

ACCESSORIES 15

MEDIGEL Z® AND MEDIPATCH GEL Z® 16


GARMENT SIZING CHART 18

KEY :


S/ Standard Line
EC/ Elegance Coolmax® Line

LIPO-PANTY® Product name
S/002 Product Reference
Waist Hight Features
● ○ ● Colors available

NEW DESIGN FOR FAT TRANSFER WITH HOLE OR BBLXL



HOLE



BBLXL



S/036

LIPOFILLING ADJUSTABLE SHOULDER

now available on demand for all Medical Z lipo Panty and Combi garment

WHO ARE WE?

For over 30 years, Medical Z[®] has strived to be at the forefront of compression garment technology.


Medical Z[®] products are hand made in USA (in Houston, TX) with the highest quality textile materials.

Thanks to its integrated R & D laboratory, Medical Z[®] has the necessary tools to control and manufacture innovative products. A requirement that has earned this family business a worldwide reputation as the front runner in quality and service specializing in burns and plastic surgery.

Medical Z[®] offers a wide range of compression orthosis (compression garments, chin straps, sternal straps) to provide patients with safety and comfort throughout the entire recovery process.



Medical Z[®] designs and manufactures products in France and the USA.



4


INTRODUCING MEDICAL Z[®]

MEDICAL Z[®] OFFERS INNOVATIVE AND SCIENTIFICALLY TESTED PRODUCTS

Our garments respect the anatomy of each patient:

The garments allow effective and constant compression even in difficult areas like the crease of the groin. They provide comfort, support, and complete range of motion.

Our high-tech fabrics consist of a fiber that provides longlasting compression, great thermoregulation, and have a high solar index of 40 - 50 SPF.



QUALITY SERVICE

Our garments are manufactured in our Houston, TX facility and can be custom made within 8 business days.

Our customer service department (English and Spanish) is available to answer any questions or comments you may have.


Toll Free : +1 800-368-7478
info.us@medicalz.com or order@medicalz.com

All materials and textiles used in Medical Z products are tested for both quality control and research in its own fabric testing laboratory based in France.

Medical Z® measures and controls the quality of its fabrics daily with state-of-the-art equipment (extensometer, perspirometer, steam equipment, and light testing equipment).

These tests include:

- Compression Force Scale
- Color Fastness to water and friction
- Dimensional stability
- Extension
- Tensile strength



MEDICAL Z® green environmental policy

In 2010, Medical Z® installed **600m² (6500 sq ft)** of photovoltaic panels on the roof of its premises in France. The production generated by solar radiation is now supplying the equivalent of **the energy consumption of Medical Z®**.

The company confirms its environmental commitment by using the **OEKO-TEX® STANDARD 100 Certified materials**.



This eco-label is the first to guarantee the **human-ecological qualities of textiles free from toxic products benefitting the body and the environment**.

The **OEKO-TEX®** system enables the rapid identification and elimination of potential sources of harmful substances for all textile products at all stages of processing.

By manufacturing its textiles with OEKO-TEX® certified raw materials, Medical Z® ensures an efficient and superior product safety.



MEDICAL BRA
EZ/001



MEDICAL BRA
POST-MASTECTOMY
EZ/005



BRA WITH SOFT FABRIC POCKETS
FOR BREAST PADS*



* Not supplied with the product.



ZBRA®
S/024

Adjustable shoulder with
mammary strap



MEDICAL ZBRA®
S/014



ZBRA®
S/036

Lipofilling adjustable shoulder



7
BREAST SURGERY



BODY HIGH BACK
EC/015

Custom made only



ZBRA®
S/021

Adjustable shoulder strap



SIZING ZBRA®
S/020



MAMMOSTAT
AC/006



MAMMARY STRAP
S/013





ALL MODELS ARE MADE WITH HOOK & EYES

*Zipper option available

BODY CONTOURING SURGERY



LIPO-PANTY®
S/001

High, calf



SHORT LIPO-PANTY®
S/004

Waist height, above knee



HOLE



SHORT LIPO-PANTY®
S/003

High, above knee



LIPO-PANTY®
S/002

Waist height, calf



BBLXL

OPTIONS NOW AVAILABLE ON DEMAND



ALL MODELS ARE MADE WITH HOOK & EYES
*Zipper option available



LIPO-PANTY®
EC/002

High, calf



LIPO-PANTY®
EC/004

High, ankle



SHORT LIPO-PANTY®
EC/003

High, above knee



SHORT LIPO-PANTY®
EC/006

Waist height, above knee



LIPO-PANTY®
EC/005

Waist height, ankle



LIPO-PANTY®
EC/001

Waist height, calf



HOLE

OPTIONS NOW AVAILABLE ON DEMAND



BBLXL

LIPO-PANTY®, PANTY-BRA® AND LIPO-BRA® postsurgical

ELEGANCE COOLMAX® LINE



ALL MODELS ARE MADE WITH HOOK & EYES

*Zipper option available

ADJUSTABLE SHOULDERS



SHORT LIPO-PANTY®

EC/034

High back, above knee



LIPO-PANTY®

EC/035

High back, calf



10

BODY CONTOURING SURGERY



COMBI-SUIT

EC/016



HOLE



ARM SLEEVES

EC/042



BBLXL


OPTIONS NOW AVAILABLE ON DEMAND

WWW.MEDICALZ.COM



PANTY-LIFT®
EC/017

Closed crotch
High, calf



ALL MODELS ARE MADE
WITH HOOK & EYES

*Zipper option available



COMBI®
EC/010

Ankle



COMBI®
EC/008

Calf ● ○



COMBI
EC/009


Above knee



HOLE

OPTIONS NOW AVAILABLE ON DEMAND

WWW.MEDICALZ.COM



BBLXL



Z-BELT® HALF TORSO FEMALE
S/006

Lateral velcro closure
Foam padding



Z-BELT® FULL TORSO FEMALE
S/008

Lateral hook closure or velcro
closure
Foam padding



ABDOMINAL BINDER FEMALE
S/022

Height: 12.5"



ABDOMINAL BINDER WITH ELASTIC FEMALE
S/022E

Height: 9"



SINGLE SIZE BINDER
S/037

Height: 9" - 12"





ALL MODELS ARE MADE WITH HOOK & EYES
*Zipper option available



BELT FULL TORSO
EC/012

Lateral hook closure



BELT
EC/033

Central closure
High back



BELT HALF TORSO
EC/011



BODY
EC/007





**ABDOMINAL BINDER
MEN**

S/022

Height: 12.5"



**ABDOMINAL BINDER
WITH ELASTIC MEN**

S/022E

Height: 9"



ABDOMINAL BINDER

S/015

Lateral velcro closure



COMBI

EC/037



ROMEO

EC/038



ROMEO VEST FOR MALE

EC/045



**Z-BELT® HALF TORSO
MALE**

S/005

Lateral velcro closure

Foam padding



**Z-BELT® FULL TORSO
MALE**

S/007

Lateral hook closure or velcro closure

Foam padding



HI-DEF FOAM FOR MEN & WOMEN



HI-DEF FOAM FOR FEMALE
S/201



HI-DEF FOAM FOR MALE
S/201



OTHER GARMENTS



VAN DE VEN
FACEMASK
AC/012

One size



VAN DE VEN
COOLING FACEMASK
AC/013

Slot for gel bag
(not supplied
with the product)
One size



EAR BAND
AC/003

Size for baby, children
and adult.
With or without foam



EAR BAND
AC/016

Velcro closure
One size



PRESSLIFT
AC/002



STERNAL STRAP
AC/001



KNEE SUPPORT
AC/004



ANKLE SUPPORT
AC/005



SLEEVE
EC/044



MEDIGEL Z[®] TECHNOLOGY

Medical Z[®] has developed a line of products that are known for their effectiveness, durability, ease in application, comfort and value.

Various studies have shown the effectiveness of polymers to prevent the formation of hypertrophic and keloid scars.


Medigel Z[®] is a combination of a triblock copolymer and an H1 Class mineral oil.

Originally devised for the treatment of scars for burn victims, Medigel Z[®] now has a wider range of products for surgeries that includes plastic surgery, aesthetics, orthopedics and podiatry.


16

MEDIGELZ[®] AND MEDIPATCH GELZ[®]


BEFORE TREATMENT



AFTER 3 MONTHS



AFTER 6 MONTHS



FINGER SLEEVES
MEDIGEL Z[®]



ELASTIC SLEEVES
MEDIGEL Z[®]



GLOVE LINED
WITH GEL Z[®]



CLEAR SHEET
MEDIGEL Z[®]

MEDIPATCH GEL Z[®] TECHNOLOGY

Medipatch Gel Z[®] has been constructed from the same technology as the Medigel Z[®], offers the advantage of self adhering to various body surfaces and has been proven superb in repetitive usage without losing its adhesive properties.

Medipatch Gel Z[®] is available as transparent or lined with fabric. It is available in various shapes and sizes to accommodate various anatomic surfaces.



MAMMOPATCH
GEL Z[®]



ABDOPATCH
GEL Z[®]



MEDIPATCH
GEL Z[®] PADS
5 X 8 CM




MEDIGEL Z® SHEETS (5 MM)	SPECIALTY PRODUCTS MEDIGEL Z®	MEDIPATCH GEL Z® SHEETS (3 MM OR 1 MM)	MEDIPATCH GEL Z® CLEAR SHEET (1 MM)	SPECIALTY PRODUCT MEDIPATCH GEL Z® (3 MM OR 1 MM) FABRIC LINED)
4 x 4 in 10 x 10 cm (Clear sheet, 2/pack)	Press-Lift (Adult or children)	2 x 3¼ in 5 x 8 cm (Fabric lined 4/pack)	2 x 3¼ in 5 x 8 cm (Fabric lined 6/pack)	Medipatch Pads with or without hole (4/pack)
10 x 10 in 20 x 20 cm (Fabric lined)	Finger Sleeves® (Adult or children, 4 sizes, 4/pack, Clear or fabric lined)	4 x 6 in 10 x 15 cm (Fabric Lined, Clear Sheet)	4 x 6 in 10 x 15 cm	Abdopatch Gel Z®
12 x 16 in 30 x 40 cm (Fabric lined)	Gloves (3 sizes)	8 x 8 in 20 x 20 cm (Fabric lined)	8 x 8 in 20 x 20 cm	Cesarpatch Gel Z®
16 x 20 in 40 x 50 cm (Fabric lined)	Elastic Sleeves (5 Sizes)	12 x 16 in 30 x 40 cm (Fabric lined)		Mammopatch Gel Z® (Left and right)
	Sternal strap (Male or female)	16 x 20 in 40 x 50 cm (Fabric lined)		

- Soft and comfortable
- Protects scar site against chafing
- Can be used in conjunction with compression garments
- Expands and conforms to patient movements
- Available with fabric lining
- Can be cut to desired size
- Hypoallergenic
- Economical: lasts 3 to 6 times longer than silicone sheets i.e. 8-12 weeks with proper care



MAIN SIZES:



- HEAD** 0a
Under the chin to the top of the head
- NECK** 0b
Middle of the neck
- THORAX CIRCUMFERENCE** 1
At axilla level
- CHEST CIRCUMFERENCE** 2
Fullest part of the nipple level
- UNDER BREAST** 3

- 4 WAIST CIRCUMFERENCE**
Smallest part of natural waist
- 5 HIP CIRCUMFERENCE**
Fullest part of hip area
- 6 THIGH CIRCUMFERENCE**
Fullest part
- 7 BICEPS CIRCUMFERENCE**
Fullest part

18

SIZING CHART

CHART



ZBRA®

S/014, S/021, S/024, S/036 (available from 30 to 38 B C D), EZ/001 (available from 32 to 42 B C D E), EZ/005 (available from 32 to 42 A B C D E)

U.S. SIZE	30	32	34	36	38	40	42	44
UNDER BREAST CIRC. (in) ③	26¾ - 28¼	28¾ - 30¼	30¾ - 32¼	32¾ - 34¼	34½ - 36¼	36½ - 38¼	38½ - 40¼	40½ - 42
CHEST CIRCUMFERENCE (in) ②								
CUP A	30¼ - 31	32¼ - 33	34¼ - 35	36¼ - 37	38¼ - 39	40¼ - 41	42¼ - 43	44¼ - 45
CUP B	31 - 32	33 - 34	35 - 36	37 - 38	39 - 40	41 - 42	43 - 44	45 - 46
CUP C	32 - 32¾	34 - 34¾	36 - 36¾	38 - 38¾	40 - 40¾	42 - 42¾	44 - 44¾	46 - 46¾
CUP D	32¾ - 33½	34¾ - 35½	36¾ - 37½	38¾ - 39½	40¾ - 41½	42¾ - 43½	44¾ - 45½	46¾ - 47½
CUP E	33½ - 34¼	35½ - 36¼	37½ - 38¼	39½ - 40¼	41½ - 42¼	43½ - 44¼	45½ - 46¼	47½ - 48¼

MAMMARY STRAP

S/013

SIZES	THORAX (in) ①
0 - 2X-Small	31½ - 33½
1 - X-Small	33½ - 35½
2 - Small	35½ - 37½
3 - Medium	37½ - 39½
4 - Large	39½ - 41¼
5 - X-Large	41¼ - 43¼

NOTE:

- For Zbra® the sizing chart match to the size **AFTER OPERATION**.
- For the other products, sizing charts match to the size **BEFORE OPERATION**.
- Custom Garments can be made
Contact us at order@medicalz.com

LIPO-PANTY, COMBI & COMBI SUIT*

S/001, S/002, S/003, S/004,

EC/001, EC/002, EC/003, EC/004, EC/005, EC/006, EC/008, EC/009, EC/010, EC/016, EC/017, EC/034, EC/035



Only for EC/016

BODY & BELT ELEGANCE COOLMAX®

EC/007, EC/011, EC/012, EC/033



SIZES	UNDER BREAST CIRC. (in) ③	WAIST CIRC. (in) ④	HIP CIRC. (in) ⑤	THIGH (in) ⑥	BICEPS CIRC. (in) ⑦
1 - 2X-Small	25½ - 29½	23¾ - 25½	32½ - 34½	17 - 18¾	9¼ - 10
2 - X-Small	29½ - 31½	25½ - 27½	35 - 36½	18¾ - 20	10 - 10¾
3 - Small	31½ - 34½	27½ - 29½	37 - 39	20 - 21½	10¾ - 11½
4 - Medium	34½ - 35½	29½ - 31½	39½ - 41¼	21 - 23½	11½ - 13
5 - Large	35½ - 40½	31½ - 36¾	41¾ - 46¾	23½ - 25½	13 - 14½
6 - X-Large	40½ - 45	37 - 42¼	47¼ - 51¼	25½ - 27	14½ - 16¼
7 - 2X-Large	45 - 51½	42½ - 47¼	51½ - 55½	27 - 28	16¼ - 17¾
8 - 3X-Large	51½ - 57¼	47½ - 51¼	56 - 59	28 - 29½	17¾ - 19¼
9 - 4X-Large	57¼ - 61¾	51½ - 55	59½ - 62½	29½ - 31	19¼ - 21¼

SIZES	HIP CIRC. (in) ⑤	WAIST CIRC. (in) ④
1 - 2X-Small	32½ - 34½	26¼ - 28¼
2 - X-Small	35 - 36½	28¾ - 30¾
3 - Small	37 - 39	31 - 33
4 - Medium	39½ - 41¼	33½ - 35¾
5 - Large	41¾ - 46¾	36 - 40
6 - X-Large	47¼ - 51¼	40½ - 44
7 - 2X-Large	51½ - 55½	44½ - 48
8 - 3X-Large	56 - 59	48½ - 52
9 - 4X-Large	59½ - 62½	52½ - 56

COMBI FOR MEN & ROMEO VEST

EC/037, EC/045

SIZES	CHEST CIRC. (in) ②	HIP CIRC. (in) ④
1 - X-Small	37½ - 39½	35½ - 37½
2 - Small	39½ - 41¼	38¾ - 39¾
3 - Medium	41¼ - 43¼	40 - 42
4 - Large	43¼ - 45¼	42½ - 44½
5 - X-Large	45¼ - 47¼	44¾ - 46½
6 - 2X-Large	47¼ - 49¼	46¾ - 49

ROMEO FOR MEN

EC/038

SIZES	CHEST CIRC. (in) ②
0 - 2X-Small	35½ - 37½
1 - X-Small	37½ - 39½
2 - Small	39½ - 41¼
3 - Medium	41¼ - 43¼
4 - Large	43¼ - 45¼
5 - X-Large	45¼ - 47¼
6 - 2X-Large	47¼ - 49¼

ABDOMINAL BINDER FOR MEN

S/015

SIZES	HIP CIRC. (in) ⑤
1 - X-Small	35½ - 37½
2 - Small	38¾ - 39¾
3 - Medium	40 - 42
4 - Large	42½ - 44½
5 - X-Large	44¾ - 46½
6 - 2X-Large	46¾ - 49

Z-BELT® STANDARD

S/005, S/006

SIZES	HIP CIRC. (in) ⑤
1 - Small	31½ - 37½
2 - Medium	37¾ - 43¼
3 - Large	43¾ - 51
4 - X-Large	51½ - 55

ARM SLEEVES

EC/042

EC/044, only biceps circumference ⑦

SIZES	BICEPS CIRC. (in) ⑦	UNDER BREAST CIRC. (in) ③
1 - X-Small	9¾ - 11½	27½ - 31½
2 - Small	11½ - 13	31½ - 35½
3 - Medium	13 - 14½	35½ - 37½
4 - Large	14½ - 16¼	37½ - 39½
5 - X-Large	16¼ - 17¾	39½ - 41¼
6 - 2X-Large	17¾ - 19¼	41¼ - 46

ABDOMINAL BINDER (MALE OR FEMALE)

S/022, S/022-E

SIZES	WAIST CIRC. (in) ④
1 - 2X-Small	25½ - 29½
2 - X-Small	29½ - 33½
3 - Small	33½ - 37½
4 - Medium	37½ - 41¼
5 - Large	41¼ - 45¼
6 - X-Large	45¼ - 49¼
7 - 2X-Large	49¼ - 53
8 - 3X-Large	53 - 57
9 - 4X-Large	57 - 61
Height : 9 or 12.5 inches	

PRESS-LIFT

AC/002

SIZES	NECK CIRC. (in) ⑥b	HEAD CIRC. (in) ⑥a
Infant	8½-9¾	16-18¾
Toddler	10¼-11½	19¼-21¾
Adult 1	11¾-13¾	22-24¾
Adult 2	14-16¾	25¼-26¾
Adult 3	16½-18½	27¼-28¼

HI-DEF FOAM

S/201

SIZES	MALE WAIST CIRC. (in) ④	FEMALE UNDER BREAST CIRC. (in) ③
X-Small	36	30
Small	39	32
Medium	41	34
Large	43	36
X-Large	45	38
XX-Large	47	40

STERNAL STRAP

AC/001

SIZES	UNDER BREAST CIRC. (in) ③
Small	32½ - 34½
Medium	35 - 36½
Large	37 - 39



* Lipo Panty and Combi: measures taken for patients between 5'2" and 5'8"



MEDICALZ®

CRÉATEUR FRANÇAIS
DE VÊTEMENTS COMPRESSIONS DEPUIS 1981

THE HIGHEST QUALITY FOR YOUR WELL-BEING

FRANCE

ZI N°1 - 14, rue Georges Cuvier - 37550 SAINT-AVERTIN
Tél. : +33 (0)2 47 71 33 33 - Fax. : +33 (0)2 47 71 33 34
Email : infofrance@medicalz.com

USA

10635 Richmond Avenue #100 - HOUSTON, TX 77042
Toll Free: +1 800-368-7478 - Phone: +1 832-767-5978
Fax: +1 832-519-9380
Email: info.us@medicalz.com
www.shopmedicalz.com

WWW.MEDICALZ.COM