

Bact-Trap™ Midi

products

Anaesthesia/Respiratory
Bacterial/Viral Filter.

This is Bact-Trap Midi, which set the standard by which every other breathing filter has to be measured.

Patient Safety and Comfort

Bact-Trap Midi offers the user a high efficiency device for breathing systems in a crystal clear housing for easy inspection. Conical fittings in accordance with ISO standard ensure good connection security.

This is Bact-Trap Midi:

- ▶ Incorporates a > 99.99% effective Bacterial/Viral Filter to prevent the spread of infectious materials into the circuits, machines and manual resuscitators.
- ▶ Isolates the patient from potential Bacterial/Viral contamination in breathing systems.
- ▶ Less frequent cleaning and decontamination of anaesthetic machines and ventilators.
- ▶ Reduces the risk of cross-infection among patients and hospital staff.
- ▶ Initial set-up with the Bact-Trap Midi Filter is simple, quick and secure.
- ▶ Different options ensure the right choice of device.
- ▶ Bact-Trap Filters should be changed every 24 hours or more frequently if resistance increase.

! Attention.
If Bact-Trap Midi is used in combination with Nebulizer or Aerosol therapy, -possible increase of resistance should be monitored.

Pharma Optimum


Different options available, see our **Design your own system**




Bact-Trap Midi, port
Code no 7110



Bact-Trap Midi port, elbow
Code no 7111



Bact-Trap Midi port, angle
Code no 7120



Product Specification

Tidal volume (ml)	100-1200
Dead space (ml) straight	49
Dead space (ml) angled	51
Weight (g)	27
Resistance (Pa) 30 l/min	110

Filtration efficiency

Bacterial	>99,99%
Fittings angled: 15 Male/15 Female - 22 Male	

Ordering Information

Bact-Trap Midi port	7110	50/400
Bact-Trap Midi port, angle	7120	50/400
Bact-Trap Midi port, elbow	7111	50/400

Disposable, for single patient use in clinically clean packaging.

INS-529v3

Distributed by:

Manufactured in Sweden by:



Pharma Systems AB
Rubanksgatan 9
741 71 Knivsta
SWEDEN

Phone: + 46 18 349500
Fax: + 46 18 349530
Email: info@pharmasys.se
Webb: www.pharmasystems.se

Bact-Trap™ **Midi**

