


FST-FB-23-180-O-O

FST-FB-23-320-O-O

Single use, sterile
colorectal biopsy forceps.

sapi med
DIAGNOSTIC
LINE




The single-use, sterile biopsy forceps with oval jaws is made of an ergonomic plastic handle with 3 rings and a 180 mm or 320 mm rigid steel shaft. It can be easily used in an office, day-surgery or even operative room setting.

TECHNICAL SPECIFICATIONS

- 3.7mm long sharp oval jaws with fenestrations.
- 3-rings handle to control the opening of the jaws.
- Compatible with all Sapi Med proctoscopes and rectoscopes.

Single use, sterile colorectal biopsy forceps.

REF FST-FB-23-180-O-O



Shaft length: 180 mm

LATEX FREE

Sterilization: ETO

Box: 10 pcs.

Manufacturer: Foresight Technology Co. LTD.


16-12, NO.18 Hongsheng Road
Jiangbei District, 400020, Chongqing, China

EC Rep.: Creo Medical S.L.U

Pol.Ind. Comarca I, C/A, N° D-22
31160 Orcoyen, Navarra, Spain
T. +34 948 29 33 24 - Fax +34 948 29 34 18

Distributed by: Sapi Med S.p.a.

REF FST-FB-23-320-O-O



Shaft length: 320 mm

LATEX FREE

Sterilization: ETO

Box: 10 pcs.

Manufacturer: Foresight Technology Co. LTD.

16-12, NO.18 Hongsheng Road
Jiangbei District, 400020, Chongqing, China

EC Rep.: Creo Medical S.L.U

Pol.Ind. Comarca I, C/A, N° D-22
31160 Orcoyen, Navarra, Spain
T. +34 948 29 33 24 - Fax +34 948 29 34 18

Distributed by: Sapi Med S.p.a.